SAFETY DATA SHEET J-647

Printed: 05/19/2015

Revision: 04/30/2015

Supersedes Revision: 04/29/2015

1. Product and Company Identification

Product Code: J-647

Product Name: Jorvet Manual & Ultrasonic Neutral pH Instrument Detergent

Company Name:Jorgensen LaboratoriesPhone Number:1450 N. VanBuren Avenue1 (800) 525-5614

1450 N. VanBuren Avenue Loveland, CO 80538

Emergency Contact: Chemtrec 1 (800) 424-9300

2. Hazards Identification

Serious Eye Damage/Eye Irritation, Category 2B

GHS Signal Word: Warning

GHS Hazard Phrases: Causes eye irritation.

GHS Precaution Phrases: Keep out of reach of children.

Read label before use. Do not get in eyes.

Wash hands thoroughly after handling.

GHS Response Phrases: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and

easy to do. Continue rinsing for 15 minutes. Get medical advice/attention.

IF ON SKIN: Wash any sensitive areas. If skin irritation occurs, get medical advice/attention.

IF SWALLOWED: Rinse mouth. Do NOT induce vomiting. Get immediate medical advice/attention. IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing. If

experiencing respiratory symptoms: Get medical advice/attention.

Medical and Transportation Emergencies contact Chemtrec 1-800-424-9300

GHS Storage and Disposal

Phrases:

Dispose of contents/container in accordance to local, state and federal

regulations.

Hazard Rating System:

Health	1
Flammability	0
Physical	0
PPE	В

HMIS:

Potential Health Effects (Acute and Chronic)

SAFETY DATA SHEET J-647

Printed: 05/19/2015

Revision: 04/30/2015

Supersedes Revision: 04/29/2015

3. Composition/Information on Ingredients

CAS # Hazardous Components (Chemical Name) Concentration

85536-14-7 Benzenesulfonic acid, 4-C10-13-sec-alkyl derivs. 10.0-20.0 %

4. First Aid Measures

Emergency and First Aid

Procedures: No data available.

In Case of Inhalation: If breathing is difficult, remove to fresh air and keep at rest in a position comfortable for breathing. If

experiencing respiratory symptoms: Get medical advice/attention.

In Case of Skin Contact: Wash thoroughly with fresh, tepid water for 15 minutes. Discard contaminated clothing and footwear or wash

before reuse. If skin irritation occurs: Get medical advice/attention.

In Case of Eye Contact: Immediately flush eyes with large amounts of fresh, tepid water for at least 15 minutes. Hold eyelids open to

ensure complete irrigation of eye and lid tissues. Tilt head to the side and irrigate the eye from the bridge of the nose to the outside of the face. Keep run-off from entering the other eye, mouth or ear. Washing eye within the first few seconds is essential to achieve maximum effectiveness. Remove contact lenses, if present

and easy to do. Continue rinsing. If eye irritation persists, get medical advice/attention.

In Case of Ingestion: Do not induce vomiting. Rinse mouth with fresh, tepid water, then immediately drink 4-8 oz. of milk or water.

Never give anything by mouth to an unconscious person. If vomiting occurs, keep airways open. Keep head

lower than hips to prevent aspiration into the lungs. Get medical advice/attention.

5. Fire Fighting Measures

Flash Pt: No data.

Explosive Limits: LEL: No data. UEL: No data.

Auto ignition Pt: No data.

Suitable Extinguishing

Media: Dry chemical, CO2, sand, earth, water spray or regular foam.

Fire Fighting Instructions: As in any fire, wear a self-contained breathing apparatus in pressure-demand, MSHA/NIOSH (approved or

equivalent), and full protective gear.

Flammable Properties and

Hazards: No data available.

6. Accidental Release Measures

Protective Precautions, Rubber or neoprene gloves, safety glasses. Wear chemical protective clothing. Do not get in eyes, on **Protective Equipment and** skin, or on clothing. Wash hands thoroughly after handling. Provide mechanical ventilation to disperse or **Emergency Procedures:** ventilate the area with fresh air.

Printed: 05/19/2015

SAFETY DATA SHEET J-647

Revision: 04/30/2015

Supersedes Revision: 04/29/2015

Steps To Be Taken in Case Absorb spill with inert material (e.g. dry sand or earth), and dispose of in accordance with applicable **Material Is Released Or** regulations.

Spilled:

7. Handling and Storage

Precautions To Be Taken For industrial or institutional use only. Avoid contact with eyes and prolonged or repeated skin contact. In Handling:

Precautions To Be Taken Keep container closed when not in use. Keep from freezing. Store away from incompatible material. **In Storing:**

8. Exposure Controls/Personal Protection

CAS # Partial Chemical Name OSHA TWA ACGIH TWA Other Limits

85536-14-7 Benzenesulfonic acid. No data. No data. No data.

4-C10-13-sec-alkyl derivs.

Respiratory Equipment: A respirator is not needed under normal and intended conditions of product use.

(Specify Type):

Eye Protection: Safety glasses

Protective Gloves: Rubber or neoprene gloves

Other Protective Clothing: Wear chemical protective clothing.

Engineering Controls Handle in accordance with good industrial hygiene and safety practice in a well ventilated area.

(Ventilation etc.): Wash hands before breaks and at the end of the workday.

9. Physical and Chemical Properties

Physical States: [] Gas [X] Liquid [] Solid Appearance and Odor: Appearance: Clear Yellow Liquid

Odor: Pleasant mint odor

Melting Point:No data.Boiling Point:No data.Autoignition Pt:No data.Flash Pt:No data.

Explosive Limits: LEL: No data. UEL: No data.

Specific Gravity (Water = 1): 1.010 - 1.030

Vapor Pressure (vs. Air or

mm Hg): No data.

Vapor Density (vs.

Air = 1): No data. Evaporation Rate: No data. Solubility in Water: 100

SAFETY DATA SHEET I-647

Printed: 05/19/2015

Revision: 04/30/2015

Supersedes Revision: 04/29/2015

pH: 6.0 – 7.0 @ 1%

Percent Volatile: N.D.

10. Stability and Reactivity

Stability: Unstable [] Stable [X]

Conditions to Avoid: Avoid handling conditions which may allow for leaks and spills of this material. Do not permit personnel to

handle this product without proper training and/or protective equipment.

Incompatibility -

Materials to Avoid: Avoid contact with strong acids, bleach bases or other.

Hazardous

Decomposition or

Byproducts: No data available.

Possibility of

Hazardous Reactions: Will occur [] Will not occur [X]

Conditions to Avoid -

Hazardous Reactions: No data available.

11. Toxicological Information

Toxicological Information: Inhalation, Eye contact, Skin contact, Ingestion Irritation or Corrosion: Eyes: Causes serious eye irritation.

Skin: No known effects.
Ingestion: No known effects.
Inhalation: No known effects.

Inna

Symptoms Related to

Toxicological Characteristics: Eyes: Redness, pain, irritation

Skin: No symptoms known.
Ingestion: No symptoms known.
Inhalation: No symptoms known.

Chronic Toxicological Effects: No known chronic toxicological effects.

Carcinogenicity: NTP? No IARC Monographs? No OSHA Regulated? No

12. Ecological Information

General Ecological Information: This product has no known ecotoxicological effects.

SAFETY DATA SHEET I-647

Printed: 05/19/2015

Revision: 04/30/2015

Supersedes Revision: 04/29/2015

13. Disposal Considerations

Waste Disposal Method: Dispose of contents/container in accordance to local, state and federal regulations.

14. Transport Information

LAND TRANSPORT (US DOT):

DOT Proper Shipping Name:

Compound Cleaner, Not Regulated

DOT Hazard Class: UN/NA Number:

AIR TRANSPORT (ICAO/IATA):

ICAO/IATA Shipping Name: Not Regulated.

15. Regulatory Information

EPA SARA (Superfund Amendments and Reauthorization Act of 1986) Lists

CAS # Hazardous Components (Chemical Name) S.302 (EHS) S.304 RQ S.313 (TRI)

85536-14-7 Benzenesulfonic acid, 4-C10-13-sec-alkyl derivs. No No No

16. Other Information

Revision Date: 04/30/2015

Additional Information About

This Product:

Company's Disclaimer: While Jorgensen believes the statements set forth herein are accurate as of the date hereof, Jorgensen makes no warranty with respect thereto and expressly disclaims all liability for reliance thereon. Such data is offered solely for your consideration, investigation, and verification.