


Covidien Needle and Syringe Selection Guide

Monoject™ Veterinary Needles & Syringes



Monoject™ Hospital Standard Needles & Syringes



KenVet™ Veterinary Needles & Syringes



The Choice is Yours

- Monoject™ Needles & Syringes, the industry leader
- KenVet™ Needles & Syringes, the economy line

Choose Your Needle

- **Veterinary Angle Point:** Designed for strength to penetrate thick or dehydrated skin
- **Hospital Point:** Designed for smooth penetration into thin well hydrated skin
- **Metal insert** for robust needle attachment
- **Epoxy needle** attachment
- **Clear or colored** hub and sheath

Choose Your Syringe

- **Luer lock or slip** tip syringe
- **Silver recycled** or white plunger

Choose Your Packaging

- **Recyclable rigid** pack
- **Peel to open** soft pack

Covidien Needle and Syringe Selection Guide

| | | Monoject™ Veterinary Needles & Syringes | Monoject™ Standard Needles & Syringes | KenVet™ Veterinary Needles & Syringes |
|----------------------------|---|--|---|--|
| Needle Sheath | Color Coded | ■ | | |
| | Clear | | ■ | ■ |
| Needle Point | Shallow needle angle point for strength to penetrate thick or dehydrated skin | ■ | | |
| | Steep needle point angle for penetration into thin well hydrated skin | | ■ | ■ |
| Needle Attachment Epoxy | Stabilizing metal insert | ■ | | |
| | Hub attachment | | ■ | ■ |
| Needle Hub | Clear | ■ | | ■ |
| | Color coded | | ■ | |
| Syringe Attachment | Luer lock | ■ | ■ | ■ |
| | Luer slip | ■ | | |
| Syringe Graduations | Heat | ■ | ■ | |
| | Printed | | | ■ |
| Finger Flanges | Low profile | ■ | ■ | |
| | Extended | | | ■ |
| Plunger Color | Silver | ■ | ■ | |
| | White | | | ■ |
| Packaging | Rigid recyclable packaging | ■ | ■ | |
| | Soft pack packaging | | | ■ |

www.covidien.com/animalhealth



COVIDIEN, COVIDIEN with logo, Covidien logo are U.S. and/or internationally registered trademarks of Covidien AG. Other Brands are trademarks of a Covidien company. ©2014 Covidien. All rights reserved.

H8202 xM 0214