

Covidien Needle and Syringe Selection Guide

Monoject™ Veterinary Needles & Syringes



Monoject™ Hospital Standard Needles & Syringes



KenVet™ Veterinary Needles & Syringes



The Choice is Yours

- Monoject™ Needles & Syringes, the industry leader
- KenVet™ Needles & Syringes, the economy line

Choose Your Needle

- **Veterinary Angle Point:** Designed for strength to penetrate thick or dehydrated skin
- **Hospital Point:** Designed for smooth penetration into thin well hydrated skin
- **Metal insert** for robust needle attachment
- **Epoxy needle** attachment
- **Clear or colored** hub and sheath

Choose Your Syringe

- Luer lock or slip tip syringe
- Silver recycled or white plunger

Choose Your Packaging

- Recyclable rigid pack
- Peel to open soft pack

Covidien Needle and Syringe Selection Guide

		Monoject™ Veterinary Needles & Syringes	Monoject™ Standard Needles & Syringes	KenVet™ Veterinary Needles & Syringes
Needle Sheath	Color Coded	<input checked="" type="checkbox"/>		
	Clear		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Needle Point	Shallow needle angle point for strength to penetrate thick or dehydrated skin	<input checked="" type="checkbox"/>		
	Steep needle point angle for penetration into thin well hydrated skin		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Needle Attachment Epoxy	Stabilizing metal insert	<input checked="" type="checkbox"/>		
	Hub attachment		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Needle Hub	Clear	<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>
	Color coded		<input checked="" type="checkbox"/>	
Syringe Attachment	Luer lock	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
	Luer slip	<input checked="" type="checkbox"/>		
Syringe Graduations	Heat	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
	Printed			<input checked="" type="checkbox"/>
Finger Flanges	Low profile	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
	Extended			<input checked="" type="checkbox"/>
Plunger Color	Silver	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
	White			<input checked="" type="checkbox"/>
Packaging	Rigid recyclable packaging	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
	Soft pack packaging			<input checked="" type="checkbox"/>

www.covidien.com/animalhealth



COVIDIEN, COVIDIEN with logo, Covidien logo are U.S. and/or internationally registered trademarks of Covidien AG. Other Brands are trademarks of a Covidien company. ©2014 Covidien. All rights reserved.

H8202 xM 0214