


Large Animal I.V. Extension Set

Large Animal 8" Extension Set with Y-Injection Site

Reorder No. EY3008

Sterile, disposable, 8 inch large bore extension set with on/off clamp, proximal female Luer lock connector, and distal rotating SPIN-LOCK™ male Luer connector with Y-injection site for use between I.V. catheter and STAT Large Animal I.V. Set.



INTERNATIONAL WIN, LIMITED

340 North Mill Road, Suite 6, Kennett Square, PA 19348-2853 USA • (610) 444-0170 • (800) 359-4946 • Fax (610) 444-0171
e-mail: info@internationalwin.com • www.internationalwin.com