

SutureVet™ Skin Stapler

Available in both 35 Wide and 35 Regular

The SutureVet™ Skin Stapler is a ready-to-use sterile disposable stapler for quick and easy skin closure.

- ☑ Ready-to-use sterile disposable stapler
- ☑ For quick and easy skin closure
- ☑ Comfortable squeeze handle design
- ☑ Good visibility allowing for accurate staple placement
- ☑ Contains 35 / 316 L stainless steel in either the Regular or Wide stapler design
- ☑ Single use/ Non-reloadable



SKU#	DESCRIPTION	PKG
029286	Staples/Clips - SutureVet Skin Stapler 35 Regular	EACH
029287	Staples/Clips - SutureVet Skin Stapler 35 Wide	EACH

VetClose® Surgical Glue

VetClose® Surgical Glue is for closure of external wounds. It is indicated for use on minor cuts, wounds, and incisions.

- ☑ High purity which gives better handling and storage
- ☑ Comes with twelve pipettes for precise administration
- ☑ Purple in its liquid state for easier visibility, but dries clear
- ☑ Strongly bonds the wound, but can flex with the wound
- ☑ Works in 5-15 seconds, depending upon the degree of moisture at the site



**For more
information on
these products
see our Surgical
Suite brochure.**



SKU#	DESCRIPTION	UOM	PKG
031477	Tissue Adhesives - SutureVet Vetclose	2mL	BOTTL



Helping to improve your clinic's profitability, without lowering your quality of care.

PVL00032-6/15