

IV Catheter

Intravenous Ethylene Tetrafluoroethylene (ETFE)



- Designed to perform in delicate vasculature, delivering maximum efficiency and patient comfort
- Smooth catheter-to-needle transition results in easy insertion
- ETFE catheter material is biocompatible, minimizing risk of irritation and infection while catheter is indwelling
- Hard plastic outer packaging protects catheter before use
- Radiopaque catheter features flashback “window” for quick visualization of blood
- Hub is color-coded for rapid identification of catheter size

ORDERING INFORMATION

SKU#	ITEM #	DESCRIPTION	PACKAGING
035037	CI+1450	14GX2	
029524	CI+1650	16GX2	
029515	CI+1832	18GX1.1/4	50/Box
029493	CI+1850	18GX2	4 Boxes/Case (200)
029516	CI+2025	20GX1	

SKU#	ITEM #	DESCRIPTION	PACKAGING
029494	CI+2032	20GX1.1/4	
029495	CI+2225	22GX1	50/Box
029497	CI+2420	24GX3/4	

CAUTION: Federal law restricts these devices to sale by or on the order of a physician.