


Foal/Calf or Multi-Use I.V. Set

Reorder No. IV8001S

Sterile, disposable, primary I.V. set with 20 drop/ml vented drip chamber (with protective fluid spike cap on vent port), 360 degree inline swivel, security tab to secure fluid spike to fluid container, two Y-injection sites, roller flow control clamp, on/off clamp, 10 foot polyurethane coil, 14 inches of tubing beyond coil, releasable strap to secure set to patient, and distal SPIN-LOCK™ male Luer connector. Set is fully compatible with the GYRO Fluid Hanger and SUPREME Ball Bearing Swivel (HS7003; IK7000).



INTERNATIONAL WIN, LIMITED

340 North Mill Road, Suite 6, Kennett Square, PA 19348-2853 USA • (610) 444-0170 • (800) 359-4946 • Fax (610) 444-0171
e-mail: info@internationalwin.com • www.internationalwin.com