

Vet Pro Infusion Pump

"They say great things come in small packages"

The Vet Pro 2000 is specifically designed to be extremely small, but very powerful for today's veterinary clinic.

User-friendly Numeric Keypad, Full Text LCD


Li-Ion Battery -15 hrs on one charge

Lightweight only 360 gm = < 1 Lb

Adjustable pressure settings

All types of IV solutions even **Blood**


Any type if IV Line can be used Performs best with Abbott products




New Piston Pumping Technology for reliability and accuracy +/- 5%

This amazing little pump is lightweight, weighing less than 1 lb and about the size of a soda. The Vet Pro 2000 can operate on A/C power while charging or remove from charging base to operate on battery alone. This will allow for greater portability making this pump the very last pump you will have to or want to buy!

settings.





Vet Pro Twin Infusion Pump

"If you think one is good, two is always better"



The Vet Twin enables you to have the ability to run dual lines at one time on the same pump!

Why tangle with multiple pumps when you can use just ONE!



Vet Twin has all the same specifications as the Vet Pro but is set up for running dual IV sets. This pump will operate off of A/C power while charging battery or on battery for approx 7.5 hrs with a charge time of only 4-6 hours. Extremely small $6 \frac{3}{4}$ " L x $4 \frac{1}{2}$ " H x 3" D and lightweight < 1 LB will take care of all clinic needs in one shot!

Jur Ek

Specifications of the Jorvet Vet Pro 2000 Infusion pump

Warranty
1 year limited warranty
Flow rate
0.1 ml/hr to 1200 ml/hr

Types of infusions Continuous, intermittent, TPN, blood

Accuracy +/- 5%

KVO Rate 0.1 ml/hr to 5 ml/hr (user adjustable)

Air Sensor Ultrasonic

Pressure settings High, Normal, Low- adjustable. Max=10PSI

Power supply 110 - 240 VAC, 50/60 Hz, 20 W

Battery Li - Ion 7.2 V, 1800 mA Rechargeable

Battery Charging Automatic when latched in A/C charger

Charge Time 4-6 hours for total charge

Alarms End Program, Low/End battery, Air inline

Door Open, Lock keypad

Mounting IV Pole mounting or Cage Clip

Dimensions $6^{3/4}$ "L x $3^{1/2}$ "W x $3^{3/8}$ "D

Weight 360 gm (11.6 oz) including mount/charging base

Catalog # J1097 Vet Pro 2000

J1097M Vet Pro MRI Compatible

J1098 Dual Infusion pump

Made in Israel

Jorgensen Laboratories Inc. 1450 N. van Buren Ave Loveland, CO 80538 1-800-525-5614 Jorvet.com