
Covidien Needle and Syringe Selection Guide



The Choice is Yours

- Monoject™ Needles & Syringes, the industry leader
- KenVet[™]Needles & Syringes, the economy line

Choose Your Needle

- Veterinary Angle Point: Designed for strength to penetrate thick or dehydrated skin
- Hospital Point:

 Designed for smooth penetration into thin well hydrated skin
- Metal insert for robust needle attachment
- Epoxy needle attachment
- Clear or colored hub and sheath

Choose Your Syringe

- Luer lock or slip tip syringe
- Silver recycled or white plunger

Choose Your Packaging

- Recyclable rigid pack
- Peel to open soft pack

Covidien Needle and Syringe Selection Guide

		Monoject [™] Veterinary Needles & Syringes	Monoject™ Standard Needles & Syringes	KenVet™ Veterinary Needles & Syringes
Needle Sheath	Color Coded			
	Clear			
Needle Point	Shallow needle angle point for strength to penetrate thick or dehydrated skin			
	Steep needle point angle for penetration into thin well hydrated skin			
Needle Attachment Epoxy	Stabilizing metal insert			
	Hub attatchment			
Needle Hub	Clear			
	Color coded			
Syringe Attachment	Luer lock			
	Luer slip			
Syringe Graduations	Heat			
	Printed			
Finger Flanges	Low profile			
	Extended			
Plunger Color	Silver			
	White			
Packaging	Rigid recyclable packaging			
	Soft pack packaging			

www.covidien.com/animalhealth



COVIDIEN, COVIDIEN with logo, Covidien logo are U.S. and/or internationally registered trademarks of Covidien AG. Other Brands are trademarks of a Covidien comapny.

©2014 Covidien. All rights reserved.

H8202 xM 0214