SPRINT BOLT Ultralight DR[™]





SPRINT BOLT Ultralight DR[™]

Speed & Flexibility Without Sacrifice

SPRINT BOLT allows you to focus on patients, not wrestle with cables or complex workflows.

You spoke we listened. The **SPRINT** BOLT Ultralight DR[™] sets up quickly, is reliable and rugged enough to handle the daily challenges of the mobile equine practice. **SPRINT** BOLT incorporates the latest 100%

wireless technologies including Wi-Fi acquisition computer, wireless Cesium panel and AutoSense.

The groundbreaking AutoSense technology automatically detects X-rays from any generator and does away with sync wires, sync boxes and sync batteries delivering 100% wireless freedom.

Powered by MUSICA™

MUSICA[™] delivers consistently high image quality across species, views and users. All of this is achieved independent of tuning files and window level adjustment.

Self-Contained


SPRINT BOLT Ultralight DR[™] self-contained digital radiography system includes acquisition unit, DR panel, wireless hardware and battery.

Self-Sustained

SPRINT BOLT can run on complete battery power or be plugged in for those endless days.

Fast

With fast acquisition time and wireless workflow options, *SPRINT* BOLT works at your pace.





Lightning Fast Setup and Workflow

Flexible

- Battery and plug in
- Remote laptop for acquisition
 and image review
- Easy transport. Multiple carrying options
- DR panel with multiple handle options

Ultralight

- 18.6 pounds self-contained, including DR panel
- 13.5 pounds horse-side acquisition laptop and wireless DR panel

Productive

Work anywhere, stow anywhere. **SPRINT** BOLT is a value priced workhorse making it a profitable asset in any practice.



Your Workflow. Your Choice.

Equine Digital Radiography Options to Fit Your Practice

Sound[™] portable digital radiography systems comprise the broadest portfolio of products in the industry.

Whether your practice is ambulatory, hospital based or mixed, you will find a DR solution that fits right in with your workflow needs.

SPRINT BOLT Ultralight DR[™]

Speed and Flexibilty

- 18.6 Pounds
- 10"x12" Cesium Wireless Panel
- Auto-sense: (No sync cables, boxes or wires)
- Battery Powered
- 4-Hours Continuous Shooting
- Multiple Carrying Options
- MacBook Pro with Retina Display Wireless Laptop



The completely wireless 10"x12" cesium sensor panel makes image acquisition easier than ever before. The removable user-specific ergonomic handle can be positioned for vertical or horizontal orientation.

YOUR IMAGE IS OUR BUSINESS.



800.268.5354

soundvet.com

