

GlucoGauge Blood Glucose Monitoring Kit

Henry Schein Animal Health is proud to offer our species specific Vet GlucoGauge Blood Glucose Meter as another high-quality, cost-effective alternative to the market's leading manufacturer brands.

Specifically calibrated for cats and dogs with an easy one button switch between species – no species coding required

- It's the only pet meter on the market with positive visual confirmation of the species being tested
- Includes species-specific results and verification during testing, in the results screen and in memory
- Only a small sample size is needed (.3uL) and get test results in 5 seconds
- Accurate and precise results: Relative difference of 3.4% for dogs and -.2% for cats compared to the leading brand relative difference of 8.4% for dogs and 9.5% for cats*

SKU#	DESCRIPTION	UOM
056198	GlucoGauge Blood Glucose Meter Kit	Each
056199	GlucoGauge Blood Glucose Strips	50 Count

*Relative differences versus the Yellow Springs International Y2300

IMPORTANT NOTE: The GlucoGauge only uses GlucoGauge Test Strips and the GlucoGauge Test Strips will not work in other meters.



KIT INCLUDES

- 1 GlucoGauge Blood Glucose Meter
- 1 GlucoGauge Owner's Manual
- 1 GlucoGauge Quick Reference Guide
- 25 GlucoGauge Blood Glucose Test Strips
- 1 GlucoGauge Control Solution
- 1 Lancing Device
- 30 Lancets
- 1 Carrying Case
- 2 CR2032 Batteries





Visit us at www.henryscheinvet.com, or contact your Henry Schein Animal Health Representative to try the Henry Schein Brand difference today!

© 2018 Henry Schein Animal Health. All Rights Reserved. No copying without permission. Not responsible for typographical errors.

