

OLYMPIC VERSA-LIFT™

One electric lift that does 3 jobs

It's a GURNEY

Light weight, easy to roll.
Turns in tight spaces.



It's a TREATMENT TABLE

Adjustable height.
Comfortable standing
or seated.



It's a LOW-BOY

Electric lift makes removing
patients from cars easier, safer.



*Watch video demonstration
www.OlympicVeterinary.com*



Saves money You save money when you purchase the Olympic Versa-Lift, because this one electric lift can be used in three different ways. It's like getting three products for the price of one. And since the Versa-Lift will be used every day, it will quickly pay for itself.

Height range The Versa-Lift can be adjusted to any height up to 39 inches for use as a treatment or exam table, or as a gurney. For low-boy curbside use, it can be lowered to near ground level.

Electric lift The electric lift is quiet and smooth – won't frighten your patients. Lifts 250 pounds, with a convenient up-down foot switch. The rechargeable battery holds a long lasting charge.

Lightweight – easy to roll When used as a gurney, the Versa-Lift is easy to roll and turn. Weighs only 68 pounds, less than half the weight of most other electric gurneys. For easy steering and tight spaces, there are two handles on either end. Large 5-inch locking casters.

Strong construction All structural joints are welded. Sturdy and safe for treating and transporting your patients.

Choice of stretcher tops Versa-Lift comes with lift-off stretcher tops held securely by 4 posts but easily lifted off. Both models have a Hard Stretcher top which is firm and non-slip for exams, treatments and transport. Standard model also includes Soft Stretcher top for transport (see next column.)

No liquid spills Liquids are confined and drained by a lip around the table top (bucket not furnished).

IV pole A 2-bag IV pole is included with all models of the Versa-Lift. IV pole fits into sockets on both sides.

Saves floor space Because one Versa-Lift takes the place of three pieces of equipment it saves floor space. Footprint: 49" x 23". Top surface: 51" x 24".

3-year Warranty Full 3-year Warranty covers *all parts and labor*, except parts subject to normal wear like Velcro straps and casters.

Lift-off stretchers – safer for staff and patients



Prevents back injuries Lift-off stretchers make it easier and safer to transfer heavy patients to cages, tables, and the floor.

Safer for patients Less traumatic, less chance of injury.



Soft Stretcher Plastic fabric is medical-grade Herculite®. Wipes clean. X-ray transparent. 2 Velcro straps.



Hard Stretcher For exams and treatments. Molded of rigid, reinforced polymer. Non-slip surface. Wipes clean.

Low-Boy features



Handle slides out

Curbside use The Versa-Lift is ideal for extracting patients from cars and vans and transporting them to the clinic. The electric lift can be raised to any level, whether the animal is on the floor or on a high seat, and the lift-off stretcher top can slide under the patient, reducing awkward handling and the chance of further injury.

Hidden handle For low-boy use, a handle slides out from under the Versa-Lift and locks in position for easier pushing over rough ground.

Hands-on care If desired, the lift can be raised while the patient is being transported so that an assistant can provide hands-on care.



Ordering and Information

Olympic Versa-Lift™. Model 50701. Standard model. Includes Hard Stretcher, Soft Stretcher, and IV pole. Electrical: 120V/60Hz.

Olympic Versa-Lift™. Model 50703. Economy model. Includes Hard Stretcher and IV pole. Electrical: 120V/60Hz.