


NEW

OLYMPIC TREATMENT TABLE™

A treatment table that is also an easy-roll gurney



Watch video demonstration
www.OlympicVeterinary.com

Multi-purpose The Olympic Treatment Table is a new combination treatment table and gurney with a stainless steel top. Older treatment tables are so heavy and hard to roll that they are practically stationary. This new table can be used for treating and prepping – then for quickly transporting patients to surgery and recovery.

Lightweight This new table rolls so easily because it is lightweight – weighs less than half as much as many older tables. Full size 5" casters roll effortlessly (then lock).

Push bar For quick steering in tight places.

Full-size stainless top 46" x 21" for large dogs.

Rechargeable lift Smooth, quiet. 250 pounds capacity.

Height adjustment Comfortable seated or standing.

2 foot switches To save time, there are up-down foot switches on both sides of the table.

Instrument tray A stainless steel instrument tray is included. The tray fits into sockets at all corners of the table and swivels. 2 trays optional. See next column.

IV pole A 2-bag IV pole is included, fits into sockets.

3-year Warranty Covers all parts and labor, except parts subject to normal wear, such as casters.

A good investment More for your money. Fully equipped treatment table with gurney.

Ordering and Information

Olympic Treatment Table™. Model 50702. Includes instrument tray and IV pole. Rechargeable lift: 120V/60Hz.

Olympic Treatment Table™ with two Instrument Trays. Model 50704. Includes IV pole. Rechargeable lift: 120V/60Hz.

Save time
One table for
treating, prepping,
and transporting

