

Start Saving Today with the Pureline® OC4000 Oxygen Concentrator.

Pureline®
by Supera Anesthesia Innovations

OC4000

Stop paying for tanks, refills, and deliveries. Stop wasting time on tank billing, scheduling, management, and maintenance. Instead, get your oxygen for free, the natural way—from the air.



Savings for you, safety for your patients.

Supera Anesthesia Innovations is proud to bring the revolutionary Pureline OC4000 oxygen concentrator to veterinarians in the U.S. for the first time. While similar machines have been in use in Europe for years, Supera has perfected the technology, added safety features, and now is making this money-saving equipment available here. The Pureline concentrator is FDA approved, UL tested, and backed by the industry's best warranty.

- Extracts medical-grade oxygen (USP93%) from the air and delivers it on demand up to 5 lpm
- Connects with industry standard oxygen DISS fittings
- Includes 10' high-pressure green oxygen hose
- Quiet operation at 40 dBA
- Light weight 30 lbs.
- Energy efficient; average 12 cents a day to operate
- Built specifically for anesthesia machine use
- FDA approved, UL tested
- Oxygen purity and flow safety alarms
- Made in the USA by craftsmen using state of the art manufacturing equipment



*Example Only
Anesthesia equipment and stand not included.*

**Return On Investment
can be in just months!**

**Visit www.superavet.com/ROI
to calculate YOUR actual savings!**



Learn more about how Pureline can benefit any
veterinary practice on our website:

www.superavet.com/pureline

Anesthesia Innovations Since 1993

SUPERA | ANESTHESIA
INNOVATIONS

Formerly LEI Medical

(503) 723-5068 or (877) 620-1500 toll-free
info@superavet.com

12330 S.E. Highway 212, Clackamas, OR 97015