

## Model M2000 **Deluxe Rebreathing Floor Stand Veterinary Anesthesia Machine**

The M2000 is a real workhorse, with excellent expansion, storage options...and our industry leading ten year warranty!





## The Supera M2000 comes with these features:


- Large heavy duty shelf for monitors or ventilators
- Patented safety pop off valve it saves patient lives!
- 0-4 Lpm oxygen flowmeter with mechanical stop control
- · Oxygen flush with safety restrictor and button shield
- CO2 absorber with easy change "One Seal" canister
- Dual O2 supply ready with check valves
- Easy to use non-rebreathing circuit common outlet
- Syringe holder (for inflating E/T cuffs)


## Breathing circuits and bags included:

Adult & Pediatric Unilimb rebreathing circuit Non-rebreathing circuit 0.5, 1, 2, 3 liter heavy duty bags



**Note:** Shown with optional Isoflurane Vaporizer





Includes our stainless steel I.V. Pole system



## SUPERA

**ANESTHESIA** INNOVATIONS

Formerly LEI Medical

(503) 723-5068 or (877) 620-1500 toll-free

info@superavet.com 12330 S.E. Highway 212 Clackamas, OR 97015