

Start Saving Today with the M2300 Wall Mount Veterinary Anesthesia Machine.

The M2300 is best suited for station work such as radiology or dental suites. It mounts on standard 16" centers and comes with our standard ten year warranty.



Note: Shown with optional Isoflurane Vaporizer



The Supera M2300 comes with these features:

- Patented safety pop off valve
- 0-4 Lpm oxygen flowmeter with mechanical stop control
- Oxygen flush with safety restrictor and button shield
- CO2 absorber with pressure gauge and a quick change "One Seal" canister
- Proven unidirectional valve design used for decades in human medicine
- Quick disconnect common gas outlet
- Syringe holder (for inflating E/T cuffs)

Breathing circuits and bags included:

- Adult Unilimb rebreathing circuit
- Pediatric Unilimb rebreathing circuit
- Non-rebreathing circuit
- 0.5, 1, 2, 3 liter bags



Includes the patented safety pop off valve

SUPERA
ANESTHESIA
INNOVATIONS

Formerly LEI Medical
(503) 723-5068 or
(877) 620-1500 toll-free
info@superavet.com
12330 S.E. Highway 212
Clackamas, OR 97015