Start Saving Today with the M2500 Wall Mount Veterinary Anesthesia Machine.

The M2500 is best suited for station work such as radiology or dental areas. Full 180 Degree movement, perfect for multiple wet & dry tables. Industry leading ten year warranty.









Includes the patented safety pop off valve

The Supera M2300 comes with these features:

- · Patented safety pop off valve
- 0-4 Lpm oxygen flowmeter with mechanical stop control
- · Oxygen flush with safety restrictor and button shield
- CO2 absorber with pressure gauge and a quick change "One Seal" canister
- Proven unidirectional valve design used for decades in human medicine
- Quick disconnect common gas outlet
- Syringe holder (for inflating E/T cuffs)

Breathing circuits and bags included:

Adult Unilimb rebreathing circuit Pediatric Unilimb rebreathing circuit Non-rebreathing circuit 0.5, 1, 2, 3 liter bags

Note: Shown with optional Isoflurane Vaporizer





INNOVATIONS

Formerly **LEI Medical**

(503) 723-5068 or (877) 620-1500 toll-free

info@superavet.com 12330 S.E. Highway 212 Clackamas, OR 97015