Model M4000 Rebreathing Table Top Veterinary Anesthesia Machine

The M4000 is our most compact anesthesia machine. It's ideal for table top or wet table locations where floor space is minimal. It comes with our Industry leading ten year warranty





The Supera M4000 comes with these features:

- Patented safety pop off valve
- 0-4 Lpm oxygen flowmeter with mechanical stop control
- Oxygen flush with safety restrictor and button shield
- CO2 absorber with pressure gauge and a quick change "One Seal" canister
- Proven unidirectional valve design used for decades in human medicine
- · Quick disconnect common gas outlet
- Syringe holder (for inflating E/T cuffs)

Note: Shown with optional Isoflurane Vaporizer





Includes the patented safety pop off valve

Includes our stainless steel I.V. Pole system

