



Sorbocryl is a sterile synthetic absorbable monofilament surgical suture, composed of poly(glycolide-co-caprolactone) (PGCL). PGCL sutures are available dyed with FD&C violet #2. PGCL suture meets all requirements established by the United States Pharmacopoeia (U.S.P.). It is comparable to Ethicon's Monocryl®<sup>1</sup> sutures.

Sorbocryl is indicated for use in soft tissue approximation and/or ligation, but not for use in cardiovascular or neurological tissue, microsurgery, or ophthalmic surgery.

This absorbable suture should not be used where extended approximation of tissue is required.

**Features:**

- Excellent holding strength
- Smooth tissue passage, excellent handling, exceptional suppleness, and good knot security
- Non-wicking material helps prevent bacterial migration along the suture line

**EQUINE SUTURE LINE** - Designed specifically for the Equine Surgeon

- **Bigger, stronger needle** options to accommodate the demands of large animal surgery
- **Large suture thread** options that are better at keeping the tougher skin of equines closed during the recovery process
- **Longer swaged suture** lengths enable surgeons to open fewer suture packs, helping to keep costs down

NEEDLE		SUTURE DESCRIPTION			
	Code	Length	Size	SKU	UOM
<p>Reverse Cutting NCP-1 36mm</p>	<p>NCP-1</p>	42" (107cm)	2/0	059038	12 ct

**Warning:** Sutures should not be used if package is opened or damaged or if the expiration date has been exceeded. Discard open unused suture. Do not resterilize; resterilization may alter the physical properties of this suture, which will result in an adverse patient reaction.

For complete product details, please reference full package insert.

**FOR VETERINARY USE ONLY**

<sup>1</sup> Trademark of Ethicon

Material	Monofilament
Color(s)	Violet
Profile	70% @ 7 days 40% @ 14 days 12% @ 21 days
Absorption Rate	91-120 days
Compare to	Monocryl® <sup>1</sup>



Helping to improve your clinic's profitability and your quality of care.