








Equine Suture Line – Designed specifically for the Equine Surgeon

- **Bigger, stronger needle** options to accommodate the demands of large animal surgery
- **Large suture thread** options that are better at keeping the tougher skin of equines closed during the recovery process
- **Longer swaged suture** lengths enable surgeons to open fewer suture packs, helping to keep costs down
- **Bright, fluorescent suture** colors make it easier for the surgeon to see the suture against dark animal skin, hair and bloody surgical sites

| | Needle | Suture Length | Size | UOM | PGA | Polyamid™ | Monosorb™ | Promend™ | Propylglo Pink™ | Sorbocryl™ |
|--|--------|---------------|------|-------|-------------------------|----------------------|---------------------|----------------------|------------------------|-----------------------|
| | | | | | MAY BE USED IN PLACE OF | | | | | |
| | | | | | Vicryl ¹ | Ethilon ¹ | PDS II ¹ | Prolene ¹ | Fluorofil ² | Monocryl ¹ |
|  <p>Reverse Cutting</p> | NCPXL | 60" (150 cm) | 3 | 12 ct | 059032 | | | | | |
| | | 42" (107 cm) | 2 | 12 ct | 059033 | 059046 | 059042 | 059049 | | |
|  <p>Reverse Cutting</p> | NCP | 42" (107 cm) | 1 | 12 ct | 059034 | 059045 | 059041 | | | |
|  <p>Reverse Cutting</p> | NCP-1 | 42" (107 cm) | 0 | 12 ct | 059035 | 059044 | 059040 | | | |
| | | | 2/0 | 12 ct | 059036 | 059043 | 059039 | | 059048 | 059038 |
|  <p>Taper Point</p> | NCT-1 | 42" (107 cm) | 2/0 | 12 ct | 059037 | | | | | |
|  <p>Straight Cutting</p> <p>NOT TO SCALE</p> | NKS | 30" (75 cm) | 2/0 | 12 ct | | 059047 | | | | |

¹ Trademark of Ethicon
² Trademark of CP Medical



Helping to improve your clinic's profitability and your quality of care.

