

Designed specifically for veterinary applications and available in a broad range of sizes, Henry Schein Brand soft pack disposable needles, syringes and combinations consistently perform at high levels for smooth tissue penetration and improved patient comfort.

SOFT PACK



Convenient soft “blister” packaging is easy-to-open, environmentally-friendly and allows hospitals to maximize storage space.



Needles

- Sharp tri-bevel point and inert lubrication for smooth tissue penetration
- ISO recognized color-coding for accurate, quick size identification
- Aluminum insert for maximum stability
- Sterilized by ethylene oxide
- Universal design works with all Luer Lock or Luer Slip syringes

Syringes:

- Available with Luer Lock and Luer Slip tips
- Catheter and eccentric tips available on select sizes
- Polypropylene barrel and elastomer piston for smooth plunger motion
- Bold barrel markings for easy and accurate dosage measurements
- Barrel designed to prevent accidental plunger removal
- Large flanges for improved control
- Tipped plunger head on 1mL size minimizes dead space
- Sterilized by ethylene oxide
- Piston is free of natural latex

KEY FEATURES



Quality, Reliability and Savings.

Needles



SKU#	BARREL	NEEDLE GAUGE	LUER	UOM
060768	—	18 G x 1.0	—	100 Count
060769	—	18 G x 1.5	—	100 Count
060762	—	20 G x 1.0	—	100 Count
060758	—	20 G x 1.5	—	100 Count
060760	—	21 G x 1.0	—	100 Count
060770	—	22 G x 3/4	—	100 Count
060759	—	22 G x 1.0	—	100 Count
060761	—	22 G x 1.5	—	100 Count
060767	—	23 G x 1.0	—	100 Count
060763	—	25 G x 5/8	—	100 Count
060764	—	25 G x 1.0	—	100 Count
060765	—	25 G x 1.5	—	100 Count
060757	—	27 G x 0.5	—	100 Count

Syringes



SKU#	BARREL	NEEDLE GAUGE	LUER	UOM
060733	3mL	—	Luer Lock	100 Count
060729	6mL	—	Luer Lock	100 Count
060731	12mL	—	Luer Lock	100 Count
060735	20mL	—	Luer Lock	50 Count
060730	30mL	—	Luer Lock	50 Count
060732	60mL	—	Luer Lock	50 Count
060734	1mL	—	Luer Slip	100 Count
060073	3mL	—	Luer Slip	100 Count
060736	6mL	—	Luer Slip	100 Count
060072	12mL	—	Luer Slip	100 Count
060739	35mL	—	Luer Slip	50 Count
060737	20mL	—	Eccentric	50 Count
060738	60mL	—	Eccentric	50 Count

Needle/Syringe Combinations



SKU#	BARREL	NEEDLE GAUGE	LUER	UOM
060728	1mL	25 G x 5/8	Luer Slip	100 Count
060722	3mL	20 G x 1	Luer Lock	100 Count
060721	3mL	20 G x 1.5	Luer Lock	100 Count
060727	3mL	21 G x 1	Luer Lock	100 Count
060067	3mL	22 G x 0.75	Luer Lock	100 Count
060720	3mL	22 G x 1	Luer Lock	100 Count
060723	3mL	22 G x 1.5	Luer Lock	100 Count
060726	3mL	23 G x 1	Luer Lock	100 Count
060725	3mL	25 G x 5/8	Luer Lock	100 Count
060071	3mL	20 G x 1.0	Luer Slip	100 Count
060066	3mL	22 G x 0.75	Luer Slip	100 Count
060068	3mL	22 G x 1.0	Luer Slip	100 Count
060069	3mL	25 G x 5/8	Luer Slip	100 Count

Quality, Reliability and Savings.